



Various tubes

Description	Ø Outside	Ø Inside	Length
Dilator 10F PEhd ORX Lg 1200	3.33 mm	1.36 mm	1'200
Dilator 11F PEhd ORX Lg 1200	3.67 mm	1.43 mm	1'200
Dilator 12F PEhd ORX Lg 1200	4.00 mm	1.50 mm	1'200
Dilator 3F PEhd ORX Lg 1200	1.00 mm	0.70 mm	1'200
Dilator 4F PEhd ORX Lg 1200	1.33 mm	0.85 mm	1'200
Dilator 5F PEhd ORX Lg 1200	1.67 mm	1.00 mm	1'200
Dilator 6F PEhd ORX Lg 1200	2.00 mm	1.07 mm	1'200
Dilator 7F PEhd ORX Lg 1200	2.33 mm	1.14 mm	1'200
Dilator 8F PEhd ORX Lg 1200	2.67 mm	1.21 mm	1'200
Dilator 9F PEhd ORX Lg 1200	3.00 mm	1.29 mm	1'200

Description	Ø Outside	Ø Inside	Length
Sheath PP 10F Lg 1180	4.33 mm	3.43 mm	1'180
Sheath PP 10F Lg 1180	4.00 mm	3.43 mm	1'180
Sheath PP 11F Lg 1180	4.67 mm	3.77 mm	1'180
Sheath PP 11F Lg 1180	4.33 mm	3.77 mm	1'180
Sheath PP 12F Lg 1180	5.00 mm	4.10 mm	1'180
Sheath PP 12F Lg 1180	4.67 mm	4.10 mm	1'180
Sheath PP 3F Lg 1180	2.00 mm	1.10 mm	1'180
Sheath PP 3F Lg 1180	1.67 mm	1.10 mm	1'180
Sheath PP 4F Lg 1200	2.33 mm	1.43 mm	1'200
Sheath PP 4F Lg 1200	2.00 mm	1.43 mm	1'200
Sheath PP 5F Lg 1200	2.67 mm	1.77 mm	1'200
Sheath PP 5F Lg 1200	2.33 mm	1.77 mm	1'200
Sheath PP 6F Lg 1200	3.00 mm	2.10 mm	1'200
Sheath PP 6F Lg 1200	2.67 mm	2.10 mm	1'200
Sheath PP 7F Lg 1180	3.33 mm	2.43 mm	1'180
Sheath PP 7F Lg 1200	3.00 mm	2.43 mm	1'200
Sheath PP 8F Lg 1180	3.67 mm	2.77 mm	1'180
Sheath PP 8F Lg 1180	3.33 mm	2.77 mm	1'180
Sheath PP 9F Lg 1180	4.00 mm	3.10 mm	1'180
Sheath PP 9F Lg 1180	3.67 mm	3.10 mm	1'180

Description	Ø Outside	Ø Inside	Length
Braided tube 6F Pebax 72D B20	2.00 mm	1.35 mm	1'250
Tube Pebax 35D 10F Lg 1200	3.33 mm	2.00 mm	1'200
Tube Pebax 35D 10F Lg 1200	3.33 mm	2.50 mm	1'200
Tube Pebax 35D 11F Lg 1200	3.67 mm	2.20 mm	1'200
Tube Pebax 35D 11F Lg 1200	3.67 mm	2.75 mm	1'200
Tube Pebax 35D 12F Lg 1200	4.00 mm	2.40 mm	1'200
Tube Pebax 35D 12F Lg 1200	4.00 mm	3.00 mm	1'200
Tube Pebax 35D 3F Lg 1200	1.00 mm	0.60 mm	1'200
Tube Pebax 35D 3F Lg 1200	1.00 mm	0.75 mm	1'200
Tube Pebax 35D 4F Lg 1200	1.33 mm	0.80 mm	1'200
Tube Pebax 35D 4F Lg 1200	1.33 mm	1.00 mm	1'200
Tube Pebax 35D 5F Lg 1200	1.67 mm	1.00 mm	1'200
Tube Pebax 35D 5F Lg 1200	1.67 mm	1.25 mm	1'200
Tube Pebax 35D 6F Lg 1200	2.00 mm	1.20 mm	1'200
Tube Pebax 35D 6F Lg 1200	2.00 mm	1.50 mm	1'200
Tube Pebax 35D 7F Lg 1200	2.33 mm	1.40 mm	1'200
Tube Pebax 35D 7F Lg 1200	2.33 mm	1.75 mm	1'200
Tube Pebax 35D 8F Lg 1200	2.67 mm	1.60 mm	1'200
Tube Pebax 35D 8F Lg 1200	2.67 mm	2.00 mm	1'200
Tube Pebax 35D 9F Lg 1200	3.00 mm	1.80 mm	1'200
Tube Pebax 35D 9F Lg 1200	3.00 mm	2.25 mm	1'200
Tube Pebax 55D 10F Lg 1200	3.33 mm	2.00 mm	1'200
Tube Pebax 55D 10F Lg 1200	3.33 mm	2.50 mm	1'200
Tube Pebax 55D 11F Lg 1200	3.67 mm	2.20 mm	1'200
Tube Pebax 55D 11F Lg 1200	3.67 mm	2.75 mm	1'200
Tube Pebax 55D 12F Lg 1200	4.00 mm	2.40 mm	1'200
Tube Pebax 55D 12F Lg 1200	4.00 mm	3.00 mm	1'200
Tube Pebax 55D 3F Lg 1200	1.00 mm	0.60 mm	1'200
Tube Pebax 55D 3F Lg 1200	1.00 mm	0.75 mm	1'200
Tube Pebax 55D 4F Lg 1200	1.33 mm	0.80 mm	1'200
Tube Pebax 55D 4F Lg 1200	1.33 mm	1.00 mm	1'200
Tube Pebax 55D 5F Lg 1200	1.67 mm	1.00 mm	1'200
Tube Pebax 55D 5F Lg 1200	1.67 mm	1.25 mm	1'200
Tube Pebax 55D 6F Lg 1200	2.00 mm	1.20 mm	1'200
Tube Pebax 55D 6F Lg 1200	2.00 mm	1.50 mm	1'200
Tube Pebax 55D 7F Lg 1200	2.33 mm	1.40 mm	1'200
Tube Pebax 55D 7F Lg 1200	2.33 mm	1.75 mm	1'200
Tube Pebax 55D 8F Lg 1200	2.67 mm	1.60 mm	1'200
Tube Pebax 55D 8F Lg 1200	2.67 mm	2.00 mm	1'200
Tube Pebax 55D 9F Lg 1200	3.00 mm	1.80 mm	1'200
Tube Pebax 55D 9F Lg 1200	3.00 mm	2.25 mm	1'200

Description	Ø Outside	Ø Inside	Length
Tube Pebax 70D 10F Lg 1200	3.33 mm	2.00 mm	1'200
Tube Pebax 70D 10F Lg 1200	3.33 mm	2.50 mm	1'200
Tube Pebax 70D 11F Lg 1200	3.67 mm	2.20 mm	1'200
Tube Pebax 70D 11F Lg 1200	3.67 mm	2.75 mm	1'200
Tube Pebax 70D 12F Lg 1200	4.00 mm	2.40 mm	1'200
Tube Pebax 70D 12F Lg 1200	4.00 mm	3.00 mm	1'200
Tube Pebax 70D 3F Lg 1200	1.00 mm	0.60 mm	1'200
Tube Pebax 70D 3F Lg 1200	1.00 mm	0.75 mm	1'200
Tube Pebax 70D 4F Lg 1200	1.33 mm	0.80 mm	1'200
Tube Pebax 70D 4F Lg 1200	1.33 mm	1.00 mm	1'200
Tube Pebax 70D 5F Lg 1200	1.67 mm	1.00 mm	1'200
Tube Pebax 70D 5F Lg 1200	1.67 mm	1.25 mm	1'200
Tube Pebax 70D 6F Lg 1200	2.00 mm	1.20 mm	1'200
Tube Pebax 70D 6F Lg 1200	2.00 mm	1.50 mm	1'200
Tube Pebax 70D 7F Lg 1200	2.33 mm	1.40 mm	1'200
Tube Pebax 70D 7F Lg 1200	2.33 mm	1.75 mm	1'200
Tube Pebax 70D 8F Lg 1200	2.67 mm	1.60 mm	1'200
Tube Pebax 70D 8F Lg 1200	2.67 mm	2.00 mm	1'200
Tube Pebax 70D 9F Lg 1200	3.00 mm	1.80 mm	1'200
Tube Pebax 70D 9F Lg 1200	3.00 mm	2.25 mm	1'200

Description	Ø Outside	Ø Inside	Length
Tube PA Grilamid L25 10F Lg 1180	3.33 mm	2.00 mm	1'180
Tube PA Grilamid L25 10F Lg 1180	3.33 mm	2.50 mm	1'180
Tube PA Grilamid L25 11F Lg 1180	3.67 mm	2.20 mm	1'180
Tube PA Grilamid L25 11F Lg 1180	3.67 mm	2.75 mm	1'180
Tube PA Grilamid L25 12F Lg 1180	4.00 mm	2.40 mm	1'180
Tube PA Grilamid L25 12F Lg 1180	4.00 mm	3.00 mm	1'180
Tube PA Grilamid L25 3F Lg 1180	1.00 mm	0.60 mm	1'180
Tube PA Grilamid L25 3F Lg 1180	1.00 mm	0.75 mm	1'180
Tube PA Grilamid L25 4F Lg 1180	1.33 mm	0.80 mm	1'180
Tube PA Grilamid L25 4F Lg 1180	1.33 mm	1.00 mm	1'180
Tube PA Grilamid L25 5F Lg 1180	1.67 mm	1.00 mm	1'180
Tube PA Grilamid L25 5F Lg 1180	1.67 mm	1.25 mm	1'180
Tube PA Grilamid L25 6F Lg 1180	2.00 mm	1.20 mm	1'180
Tube PA Grilamid L25 6F Lg 1180	2.00 mm	1.50 mm	1'180
Tube PA Grilamid L25 7F Lg 1180	2.33 mm	1.40 mm	1'180
Tube PA Grilamid L25 7F Lg 1180	2.33 mm	1.75 mm	1'180
Tube PA Grilamid L25 8F Lg 1180	2.67 mm	1.60 mm	1'180
Tube PA Grilamid L25 8F Lg 1180	2.67 mm	2.00 mm	1'180
Tube PA Grilamid L25 9F Lg 1180	3.00 mm	1.80 mm	1'180
Tube PA Grilamid L25 9F Lg 1180	3.00 mm	2.25 mm	1'180